



## IHC SELECT® BIOTINYLATED SECONDARY ANTIBODY GOAT ANTI-RABBIT IgG

---

<b>Catalog Number:</b>	21537
<b>Volume:</b>	15mL
<b>Format:</b>	Goat Anti-Rabbit IgG in PBS containing carrier protein, with 0.2% Tween 20 and 0.1% ProClin 300 as preservative.
<b>Application:</b>	<p>Prediluted for use in Immunohistochemical staining procedures. For example:</p> <ul style="list-style-type: none"><li>○ Blocking Reagent</li><li>○ Primary Antibody (Rabbit IgG) or Negative Control Reagent</li><li>○ <b>Biotinylated Secondary Antibody (Gt X Rb IgG)</b></li><li>○ Streptavidin-Enzyme Conjugate</li><li>○ Chromogen</li><li>○ Counterstain</li></ul> <p>For best results use this reagent in conjunction with other IHC Select™ detection reagents. A complete list of reagents is available at <a href="http://www.chemicon.com">www.chemicon.com</a>.</p>
<b>Storage / Handling:</b>	When stored at 2-8°C, this reagent is stable up to the expiration date printed on the label. Do not freeze or expose to elevated temperatures.

General Purpose Reagent