

77093 CCDA Selective Supplement (Charcoal Cefoperazone Desoxycholate Agar Selective Supplement)

An antimicrobial supplement recommended for the selective isolation of *Campylobacter* or *Arcobacter* species.

Composition

(per vial, sufficient for 500 ml medium):

Cefoperazone 16 mg

Amphotericin B 5 mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, cooled (45-50°C) Blood Free *Campylobacter* Broth (Cat. No. 59751). Mix well and dispense as desired.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

