



QUALITY CONTROL RANGES
MILLIPLEX[®] MAP Rat Vascular Injury Magnetic Bead Panel 2
Kit Catalog# RV2MAG-26K
Control Catalog # RV2-6026 Lot# RV2-104 and RV2-204

Analyte	QC Level	Expected Range	Units	Analyte	QC Level	Expected Range	Units
sICAM-1	Control 1	0.37 - 0.77	ng/mL	Adiponectin	Control 1	2.4 - 5.0	ng/mL
	Control 2	3.2 - 6.6	ng/mL		Control 2	20 - 41	ng/mL
sE-Selectin	Control 1	0.48 - 1.0	ng/mL				
	Control 2	4.1 - 8.6	ng/mL				

Note: The Quality Control Ranges are generated with overnight assay format.