

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

# **ProductInformation**

#### Sodium diatrizoate dihydrate

Product Number **\$4506** Store at Room Temperature

#### **Product Description**

Molecular Formula: C<sub>11</sub>H<sub>8</sub>I<sub>3</sub>N<sub>2</sub>O<sub>4</sub>Na•2H<sub>2</sub>O

Molecular Weight: 671.93 CAS Number: 737-31-5

Melting Point: 261-262 °C (with decomposition)<sup>1</sup> Synonyms: 3,5-bis(acetylamino)-2,4,6-triiodobenzoic acid sodium salt; 3,5-diacetamido-2,4,6-triiodobenzoic

acid sodium salt; sodium 3,5-diacetamido-

2,4,6-trijodobenzoate; diatrizoate sodium; Hypaque®

sodium '

Sodium diatrizoate is used as a density gradient reagent in blood cell separation. The use of sodium diatrizoate in a discontinuous gradient protocol for neutrophil separation in whole blood has been described. Sodium diatrizoate has also been utilized as a radiological contrast agent. Several studies have examined the effects of diatrizoate on mammalian cell viability, such as on cultured human vascular smooth muscle cells.

Diatrizoate has been shown to cause partial redistribution of tight-junction-associated membrane proteins into a cytoplasmic compartment in cultured monolayer MDCK cells. In cultured human umbilical vein endothelial cells, diatrizoate has been demonstrated to lead to concentration-dependent cell death and an opening of tight junctions. The interaction of diatrizoate with cyclophosphamide and carmustine in a rat bone marrow cell model has been investigated.

An investigation of the occurrence of diatrizoate in the aquatic environment has been published. <sup>7</sup>

### **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

# **Preparation Instructions**

This product is soluble in water (350 mg/ml), yielding a clear to slightly hazy, colorless to faint yellow solution. The pH of a 50% solution is 7.0 - 7.5.

Solutions of this product may be sterilized by autoclaving.

# Storage/Stability

Solutions should be protected from light.

#### References

- 1. The Merck Index, 12th ed., Entry# 3040.
- Schinella, M., Improved neutrophil separation from anti-D-treated Rh(D)-positive whole blood by a discontinuous gradient method. Am. J. Clin. Pathol., 96(3), 391-393 (1991).
- 3. Wang, Y. X., et al., The effect of radiographic contrast media on human vascular smooth muscle cells. Br. J. Radiol., **71(844)**, 376-380 (1998).
- Schick, C. S., and Haller, C., Comparative cytotoxicity of ionic and non-ionic radiocontrast agents on MDCK cell monolayers in vitro. Nephrol. Dial. Transplant., 14(2), 342-347 (1999).
- Fauser, C., et al., Differential effects of radiocontrast agents on human umbilical vein endothelial cells: cytotoxicity and modulators of thrombogenicity. Eur. J. Med. Res., 6(11), 465-472 (2001).
- 6. Parvez, Z., et al., Cytogenetic interactions of contrast media and antineoplastic drugs. Invest. Radiol., **23(5)**, 389-393 (1988).
- 7. Putschew, A., et al., Occurrence of triiodinated X-ray contrast agents in the aquatic environment. Sci. Total Environ., **255(1-3)**, 129-134 (2000).

Hypaque is a registered trademark of Sterling-Winthrop, Inc..

GCY/RXR 2/07

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.