

## QUALITY CONTROL RANGES



### MILLIPLEX<sup>®</sup> MAP Human MMP Magnetic Bead Panel 2

Control Catalog # HMMP2-6055-2  
Lot# MMP2-114 and MMP2-214

<u>Hormone</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
HMMP-1	Control I	268 - 556	pg/mL
	Control II	1324 - 2750	pg/mL
HMMP-2	Control I	686 - 1424	pg/mL
	Control II	3102 - 6442	pg/mL
HMMP-7	Control I	5887 - 12226	pg/mL
	Control II	28854 - 59927	pg/mL
HMMP-9	Control I	134 - 278	pg/mL
	Control II	785 - 1631	pg/mL
HMMP-10	Control I	271 - 563	pg/mL
	Control II	1378 - 2862	pg/mL

HMMP2-6055-2: MMP2-114.214 (14 June 2023)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • [emdmillipore.com](http://emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA