

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

METHYLPARABEN SIGMA REFERENCE STANDARD

Product No. M8911

Lot No: 111K8800 CAS No: 99-76-3

Storage Temperature: room temperature

Synonyms: Benzoic acid, 4-hydroxy-, methylester.

Methyl *p*-hydroxybenzoate.

Formula: C₈H₈O₃

Formula Weight: 152.15

Product Description

This product meets USP specifications (Current through USP 24, 4th Supplement), and is traceable to USP Reference Standard lot J


Keep tightly closed at room temperature. Dry a portion over silica gel for 5 hours before using. Use promptly. Discard unused material.

TEST	USP 24 SPECIFICATIONS	Result Sigma R.S. Lot 111K8800
Infrared absorption spectrum KBr dispersion (197M)	compares to standard	conforms to USP RS lot J
Identification B	the principal spot obtained in the chromatogram of $Test\ Solution\ B$ corresponds in size and R_f to that of the principal spot from $Standard\ solution\ B$	passes
Identification C	passes test	passes
Color of solution	passes test	passes
Acidity	passes test	passes
Residue on ignition (281)	not more than 0.1%	passes
Chromatographic purity (621)	no impurity greater than 0.5%	passes
Melting range (741)	between 125° and 128°	126° - 127°
Organic volatile impurities Method IV (467)	meets requirements	passes
Assay 〈541〉	not less than 99.0% and not more than 100.5%	100.1%

See reverse side for IR spectrum of Sigma Reference Standard.

METHYLPARABEN SIGMA REFERENCE STANDARD

Product No. **M8911** Lot 111K8800



11/2001