



**MILLIPLEX® MAP Mouse Kidney Injury
Magnetic Bead Panel 1
Kit Catalog # MKI1MAG-94K
Control Catalog # MKI1-6094-1
Lot# MKI1-108 and MKI1-208**

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
b-2-Microglobulin	Control 1	310 - 644	pg/mL
	Control 2	3184 - 6612	pg/mL
Renin	Control 1	312 - 648	pg/mL
	Control 2	3216 - 6679	pg/mL
KIM-1	Control 1	156 - 324	pg/mL
	Control 2	1612 - 3347	pg/mL
TIMP-1	Control 1	140 - 291	pg/mL
	Control 2	1590 - 3303	pg/mL
IP-10	Control 1	33 - 70	pg/mL
	Control 2	305 - 634	pg/mL
VEGF	Control 1	6.6 - 14	pg/mL
	Control 2	67 - 139	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control Values in culture media are not tested.