


DL- α -Hydroxystearic acid

Mass: 300.2664 Formula: C₁₈H₃₆O₃

mode: – collisionE: 0V

METLIN ID: 35432 SIGMA Prod: H9631

http://metlin.scripps.edu/metabo_info.php?molid=35432




DL- α -Hydroxystearic acid

Mass: 300.2664 Formula: C₁₈H₃₆O₃

mode: - collisionE: 10V

METLIN ID: 35432 SIGMA Prod: H9631

http://metlin.scripps.edu/metabo_info.php?molid=35432




DL- α -Hydroxystearic acid

Mass: 300.2664 Formula: C₁₈H₃₆O₃

mode: - collisionE: 20V

METLIN ID: 35432 SIGMA Prod: H9631

http://metlin.scripps.edu/metabo_info.php?molid=35432



DL- α -Hydroxystearic acid

Mass: 300.2664 Formula: C₁₈H₃₆O₃

mode: - collisionE: 40V

METLIN ID: 35432 SIGMA Prod: H9631

http://metlin.scripps.edu/metabo_info.php?molid=35432

