

QUALITY CONTROL RANGES



MOUSE SERUM AMYLOID A-3 ELISA KIT CATALOG # EZMSAA3-12K

STANDARD, LYOPHILIZED CATALOG # E8012-K

STANDARD	MOUSE SAA-3 VALUE
LOT NUMBER: 812-10K	5.7 µg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.25 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6012-K

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	MSAA3-110	0.12 - 0.25	µg/mL
QC 2	MSAA3-210	1.1 - 2.4	µg/mL

Reconstitute each vial with **0.25 mL** distilled or deionized water.

E6012-K: MSAA3-110.210 (13 June 2023)
EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)

www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.
© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA