



## QUALITY CONTROL RANGES

**MILLIPLEX<sup>®</sup> MAP Human Neurological Disorders Magnetic Bead Panel 3**  
**Kit Catalog # HND3MAG-39K**  
**Control Catalog # HND3-6039-3**  
**Lot# HND3-107 and HND3-207**

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
SOD1	Control I	0.30 - 0.63	ng/mL
	Control II	3.0 - 6.1	ng/mL
Angiotensinogen	Control I	10 - 21	ng/mL
	Control II	101 - 209	ng/mL
Fetuin A	Control I	32 - 67	ng/mL
	Control II	302 - 627	ng/mL
SOD2	Control I	0.24 - 0.49	ng/mL
	Control II	2.4 - 5.0	ng/mL
Contactin-1	Control I	0.62 - 1.3	ng/mL
	Control II	6.0 - 12	ng/mL
Kallikrein-6	Control I	0.29 - 0.61	ng/mL
	Control II	2.8 - 5.8	ng/mL
Osteopontin	Control I	1.1 - 2.3	ng/mL
	Control II	10 - 21	ng/mL

**Note: The Quality Control Ranges are generated with overnight assay format.  
Quality Control values in culture media are not tested.**

HND3-6039-3: HND3-107.207 (17 July 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • [emdmillipore.com](http://emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA