



Eppendorf – 5415 D and 5415 R

M8140 and Z604062 (supplied with rotor F 45-24-11), **Z605212** (supplied without rotor)

*Rotor F 45-24-11 includes polypropylene lid 5425 701 004 and rotor FA 45-24-11 includes aerosol tight aluminium lid ** Rotor F-45-36-8 supplied with the polypropylene lid 5425 733 003

ml	Tubes per rotor	RCF	Rotor Cat. No. / model	Swing out or fixed	Suspension lid Cat. No. / model	Adapters Cat. No. / model
1.5-2.0	24	16,110	* R7018 / F 45-24-11 Z666408 / FA 45-24-11	F	Polypropylene Z606693 / 5425 701 004* OR Aluminium 5425 010 000	None
0.5-0.6	24	16,110	* R7018 / F 45-24-11 Z666408 / FA 45-24-11	F	Polypropylene Z606693 / 5425 701 004* OR Aluminium 5425 010 000	Z366080 / 5425 716 001
0.4	24	16,110	* R7018 / F 45-24-11 Z666408 / FA 45-24-11	F	Polypropylene Z606693 / 5425 701 004* OR Aluminium 5425 010 000	Z366099 / 5425 717 008
0.2	24	16,110	* R7018 / F 45-24-11 Z666408 / FA 45-24-11	F	Polypropylene Z606693 / 5425 701 004* OR Aluminium 5425 010 000	A4190 / 5425 715 005
0.2 PCR tubes	36	16,110	R9272 / F-45-36-8**	F	Polypropylene Z606685 / 5425 733 003	A3566 / 5425 723 008