

THE DOZN™ SCALE

Based on the 12 Principles of Green Chemistry*, DOZN helps researchers, scientists, and manufacturers increase performance and efficiency while reducing human and environmental impact.

*Paul T. Anastas and John C. Warner, 1991.



1-Aminobenzotriazole (A3940)


	12 Principles of Green Chemistry	Percentage of Improvement	Results
Resource Used	Atom Economy	67%	Increased yield. Used less raw materials
	Waste Prevention	68%	Reduced waste by decreasing solvent usage
	Reduce Derivatives	N/A	
	Renewable Feedstocks Use	67%	Reduced auxiliary chemicals
	Real-Time Pollution Prevention	N/A	
	Catalyst	N/A	
Human & Environmental Hazards Reduction	Energy Efficiency Design	43%	Reduced chemical processing
	Less Hazardous Chemical Synthesis	36%	Eliminated hazardous hydrogenation procedure
	Safer Chemical Design	N/A	
	Safer Solvents and Auxiliaries	70%	Reduced solvent usage
	Design for Degradation	N/A	
	Inherently Safer Chemical for Accident Prevention	72%	Eliminated highly reactive chemicals

TOTAL PERCENT IMPROVEMENT

51%

AGGREGATE SCORE

0= Most Desirable



The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

© 2020 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M and DOZN are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. 2020 - 32017