


## Product Information

### 93662 Fluorescent Red 700

Fluorescent red 700 is a new fluorescent dye for the near infrared area. Extremely long wavelength fluorescence is especially well suited for applications where sample autofluorescence might be critical or penetration of tissue or other matrix is required.

#### Product Description

MW	668.86
Molar absorbance (determined in ethanol)	96.000 l · mol <sup>-1</sup> · cm <sup>-1</sup>
Abs. Max	702 nm (Ethanol)
Emission Max.	723 nm (Ethanol)
Quantity	1 mg



#### Spectrum

