

01318 *Enterococcus faecium* Selective Supplement

For selective isolation and detection of *Enterococcus faecium*.

Composition

(per vial, sufficient for 500 ml medium)

Aztreonam 37.5mg

Cephalexin 25mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, molten *Enterococcus faecium* ChromoSelect Agar Base (Cat. No 90919). Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

