


ITEM NO	DESCRIPTION	QTY	UOM
1	200L STANDARD BAG	1	EA
2	C-Flex Clear 374, 1/2" ID x 3/4" OD x 36"L	5.0	EA
3	1/2" MPX Coupling Body	5.0	EA
4	1/2" MPX Sealing Plug	5.0	EA
5	1/2" Port Plate	5.0	EA
6	Clamp-pinch 3/4" (11/16"OD to 3/4"OD tube)	5.0	EA
7	Cable Tie,7.5 White	10.0	EA
8	Silicone Pinch Clamp Covers .225 x 1.09 (Large)	10.0	EA



EMD MILLIPORE CORPORATION
BURLINGTON, MASSACHUSETTS 01803



TITLE:
MOBIUS® ESSENTIAL
ASSEMBLY, 200L, BAG

SIZE B	DWG NO ESA0200L41DP5	VER 1.0
------------------	--------------------------------	-------------------

SCALE NONE	AUTHORING SITE MM- CONTRACT MANUFACTURING	SECONDARY APPROVAL SITE N/A	DOCUMENT SECURITY N/A	SHEET 1 OF 1
---------------	---	-----------------------------------	-----------------------------	-----------------

NOTICE: This document is proprietary and its contents are the exclusive property of Merck KGaA, Darmstadt, Germany and its affiliates. This document may not be reproduced in any form whatsoever, without prior written permission from Quality Assurance Department of the life science business of Merck KGaA Darmstadt, Germany.