

Product Information

Anti-Goat IgG (Whole Molecule)–Peroxidase Antibody Produced in Rabbit

Affinity isolated antibody, Buffered aqueous solution

A4174

Product Description

Anti-Goat IgG (whole molecule) is developed in rabbit using purified goat IgG as the immunogen. Affinity isolated antigen specific antibody is obtained from rabbit anti-goat IgG antiserum by immunospecific purification which removes essentially all rabbit serum proteins, including immunoglobulins, which do not specifically bind to the goat IgG. Rabbit anti-goat IgG is then conjugated to Horseradish Peroxidase, Type VI (P8375) by a modification of the periodate method of Wilson and Nakane.¹

Specificity of Peroxidase Conjugated Anti-Goat IgG is determined by immunoelectrophoresis (IEP) versus normal goat serum and goat IgG, prior to conjugation. Cross reactivity of the antibody-conjugate is determined by ELISA. The conjugate shows no reactivity with human serum proteins.

Identity and purity of the antibody is established by immunoelectrophoresis, prior to conjugation. Electrophoresis of the antibody preparation followed by diffusion versus anti-rabbit IgG and anti-rabbit whole serum results in single arcs of precipitation.

Reagents

The conjugate is provided as a solution in Solution in 0.01 M phosphate buffered saline, pH 7.4, containing 1% bovine serum albumin and 0.05% MIT.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses.

Storage/Stability

For continuous use, store at 2-8 °C. For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Product Profile

We are now reporting lot specific information as a titer by direct ELISA (minimum 1:12,000) rather than as a working dilution. Titer is defined as the dilution of conjugate sufficient to give a change in absorbance of 1.0 at 450 nm after 30 minutes of substrate conversion at 25 °C. Microtiter plates are coated with purified goat IgG at a concentration of 5 µg/mL in 0.05 M carbonate-bicarbonate buffer, pH 9.6 (Carbonate-Bicarbonate Buffer Capsules are available as C3041).

Substrate

o-Phenylenediamine Dihydrochloride (OPD, P8287), 0.4 mg/mL in 0.05 M phosphate-citrate buffer, pH 5.0 containing 0.03% sodium perborate (Phosphate-Citrate Buffer Capsules with sodium perborate are available as P4922).

Working dilution should be determined by titration assay. Due to product improvement and changes in the assay procedure, we now list a lot specific titer by direct ELISA for this product. Due to differences in assay systems, this titer may not reflect the user's actual working dilution.

References

1. Wilson, M. and Nakane, P., Immunofluorescence and Related Staining Techniques, (Elsevier-North Holland Biomedical Press, Amsterdam), p. 215, (1978).
2. Voller, A., et al., Bull. World Health Organ., 53, 55 (1976).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2020-2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
A4174pis Rev 11/24

**MILLIPORE
SIGMA**