



## QUALITY CONTROL RANGES

**MILLIPLEX<sup>®</sup> MAP Human Complement Magnetic Bead Panel 1**  
**Kit Catalog# HCMP1MAG-19K**  
**Control Catalog # HCMP1-6019-1**

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Complement C2	Control 1	HCM1-107	0.55 – 1.1	ng/mL
	Control 2	HCM1-207	2.2 – 4.6	ng/mL
Complement C4b	Control 1	HCM1-107	11 – 22	ng/mL
	Control 2	HCM1-207	45 – 93	ng/mL
Complement C5	Control 1	HCM1-107	23 – 47	ng/mL
	Control 2	HCM1-207	90 – 186	ng/mL
Complement C5a	Control 1	HCM1-107	32 – 67	pg/mL
	Control 2	HCM1-207	123 – 256	pg/mL
Complement Factor D	Control 1	HCM1-107	0.52 – 1.1	ng/mL
	Control 2	HCM1-207	2.0 – 4.2	ng/mL
MBL	Control 1	HCM1-107	1.1 – 2.3	ng/mL
	Control 2	HCM1-207	4.3 – 8.9	ng/mL
Complement Factor I	Control 1	HCM1-107	6.0 – 12	ng/mL
	Control 2	HCM1-207	24 – 49	ng/mL

HCMP1-6019-1: HCM1-107.207 (August 17, 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA