

1.09502.0001

Supelco

MQuant®

100 Tests

pH-indicator strips (non-bleeding)


pH 2.0 - 9.0

Special indicator



for pH measurement in turbid solutions (suspensions)

In turbid solutions, on surfaces or pastes, evaluate the color through the back of the strip (for further information see the website).



7.69735.9612 Version 06/2020 - XXXXXX

Merk KGaA, 64271 Darmstadt, Germany
Tel. +49(0)6151 72-2440
EMD Millipore Corporation, 400 Summit Drive
Burlington MA 01803, USA, Tel. +1-978-715-4321
Sigma-Aldrich Canada Co. or Millipore (Canada) Ltd.
2149 Winston Park, Dr. Oakville, Ontario, L6H 6J8
Phone: +1 800-565-1400

Tiras indicadoras pH 2,0 - 9,0 (no destiñen)
para medir el pH en soluciones turbias (suspensiones)
Para las soluciones turbias, las superficies y las pastas, evaluar
el color a través del dorso de la tira (para más información,
ver sitio web).

Bandelettes indicatrices pH 2,0 - 9,0 (ne déteignant pas)
pour la mesure du pH dans les solutions troubles (suspensions)
Pour les solutions troubles, les surfaces et les pâtes, évaluer
la couleur au travers du verso de la bandelette (pour de plus
amples informations, cf. site web).

Indikatorstäbchen pH 2,0 - 9,0 (nicht blutend)
zur pH-Messung in trübten Lösungen (Suspensionen)
Bei trübten Lösungen, Oberflächen oder Pasten Farbe durch
die Rückseite des Stäbchens hindurch beurteilen (Näheres s.
Webseite).



Made in Germany

pH 2.0 - 9.0