

1.10648

Supelco



sigmaaldrich.com

100 Tests

MQuant®

Carbonate Hardness Test

Carbonathärte-Test



$1^{\circ}\text{d} \triangleq 10 \text{ mg/l CaO} \triangleq 1.25^{\circ}\text{e} \triangleq 1.78^{\circ}\text{f} \triangleq 17.8 \text{ mg/l CaCO}_3$
 $\triangleq 0.18 \text{ mmol/l CaCO}_3 (\text{Ca}) \triangleq 0.36 \text{ mmol/l HCO}_3^-$
 $\triangleq 21.8 \text{ mg/l HCO}_3^-$



For IfU please scan code with My M Safety App or go to sigmaaldrich.com



©2025

Store cool and dry • (15 - 25 °C) • Kühl und trocken lagern

Lot

min. shelf life (YYYY/MM/DD)