

QUALITY CONTROL RANGES

Mouse FGF-23 ELISA CATALOG # EZMFGF23-43K

STANDARD, LYOPHILIZED CATALOG # E8043-K

STANDARD	Mouse FGF-23 VALUE
LOT NUMBER: 843 – 2K	80 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.25mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6043-K

Reconstitute each vial with **0.25mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	MFGF23 – 102	0.48 – 1.0	ng/mL
QC 2	MFGF23 – 202	5.3 – 11	ng/mL