



QUALITY CONTROL RANGES

RAT / MOUSE GHRELIN (ACTIVE) ELISA KIT CATALOG # EZRGRA-90K

STANDARD, LYOPHILIZED CATALOG # E8090-K

Reconstitute standard with **2.0mL** distilled or deionized water.

| | |
|-----------------------------------|------------------------------------|
| STANDARD LOT NUMBER: RGA90-22K | ACTIVE GHRELIN VALUE 1980 pg/mL |
|-----------------------------------|------------------------------------|

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6090-K

Reconstitute each vial with **0.5mL** distilled or deionized water.

| CONTROL | LOT NUMBER | EXPECTED RANGE | UNITS |
|---------|------------|----------------|-------|
| QC 1 | ARGR-122 | 51 – 111 | pg/mL |
| QC 2 | ARGR-222 | 280 – 583 | pg/mL |

E6090-K: ARGR-122.222 (22 October 2024)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA