



**QUALITY CONTROL RANGES
HUMAN ADIPONECTIN RIA KIT
CATALOG # HADP-61HK**

**STANDARD, LYOPHILIZED
CATALOG# 8061-K**

STANDARD	HUMAN ADIPONECTIN VALUE
LOT NUMBER: 15HA61	200 ng/mL

The exact value of lyophilized standard is shown above.
Reconstitute with **1 mL** distilled or deionized water.

**QUALITY CONTROLS 1 AND 2, LYOPHILIZED
CATALOG #6361-K**

Reconstitute each vial with **1 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	HAD-110	6.5 – 13	ng/mL
QC 2	HAD-210	46 – 95	ng/mL

6361-K: HAD-110.210 (14 July 2021)
EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)
www.emdmillipore.com