

Specification Sheet

Product Name:	PIGMENT FREE, FILLER FREE SLEEVE STOPPER SEPTA
Product Numbers:	Z566136, Z566144, Z566152, Z566160, Z566179
Color:	AMBER
Application:	SLEEVE STOPPER
Polymer:	NATURAL RUBBER
Filler:	NONE
Cure System:	CONVENTIONAL
Formulation:	PROPRIETARY

NOMINAL PROPERTIES

HARDNESS: (SHORE A – ASTM D2240)	28 +/- 5
TENSILE STRENGTH: (PSI – ASTM D-412)	2235
MODULUS @ 100%: (PSI - ASTM D-412)	68
ELONGATION: (PERCENT - ASTM D-412)	684
SPECIFIC GRAVITY:	0.94

THIS DATA IS PRESENTED FOR INFORMATION PURPOSES ONLY.
NO WARRANTY OR RECOMMENDATION FOR ANY PARTICULAR END USE IS INTENDED OR IMPLIED.
Revised 1-25-08