



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sia.com  
sigma-aldrich.com

## Product Information

### ANTI-DOG WHOLE SERUM

Developed in Rabbit  
Delipidized, Whole Antiserum

Product No. D 4908

#### Product Description

Anti-Dog Whole Serum was developed in rabbit using pooled normal dog serum as the immunogen. The antiserum has been treated to remove lipoproteins.

Strong reactivity with normal dog serum has been determined by immunoelectrophoresis (IEP). This antiserum has not been assayed for interspecies cross reactivity.

#### Reagents

Rabbit Anti-Dog Serum is provided as a liquid containing 0.1% sodium azide as preservative.

#### Precautions and Disclaimer

Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

#### Storage and Stability

Store at 2-8 °C for up to one month.

For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

PCS/RYC 05/04

Sigma brand products are sold through Sigma-Aldrich, Inc.  
Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.  
Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.  
Please see reverse side of the invoice or packing slip.