



QUALITY CONTROL RANGES

SARS-COV2 ANTI-S1 IgG ELISA CATALOG # EZS1G-130K

STANDARD, LYOPHILIZED CATALOG # E8130-K

| | |
|---------------------|--------------------------------|
| STANDARD | SARS-COV2 ANTI-S1 IgG VALUE |
| LOT NUMBER: 8130-2K | 84 ng/mL |

The exact value of lyophilized standard is shown above.

Reconstitute with **1.0mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6130-K

Reconstitute each vial with **1.0mL** distilled or deionized water.

| CONTROL | LOT NUMBER | EXPECTED RANGE | UNITS |
|---------|------------|----------------|-------|
| QC 1 | 6130-102 | 3.2 - 6.6 | ng/mL |
| QC 2 | 6130-201 | 13 - 26 | ng/mL |

E6130-K: 6130-102.202 (17 February 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA