sigma-aldrich.com



3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

69244 Mant-GDP triethylammonium salt solution

Molecular Formula: Molecular Weight: Form: Concentration: pH: Spectroscopic Properties:

Package quantity:

Unit Definition:

 $\begin{array}{l} C_{18}H_{22}N_6O_{12}P_3 \mbox{ (free acid)} \\ 576.35 \mbox{ g/mol (free acid)} \\ \mbox{ clear aqueous solution} \\ 10 \mbox{ mM} \\ 7.5 \mbox{ +/-}0.5 \\ \lambda_{max} \mbox{ 252/355 nm}; \mbox{ } \lambda_{exc} \mbox{ 355 nm}; \mbox{ } \lambda_{em} \mbox{ 448 nm}; \\ \mbox{ ϵ 22.6/5.7 L mmol^{-1} cm^{-1} (Tris-HCl pH \mbox{ 7.5}) \\ 150 \mbox{ Units} \\ 1 \mbox{ unit = 1 ul of a 10 mM solution} \end{array}$



Structural formula of Mant-GDP

Applications:

```
Conformational dynamic: DnaB/C-proteins<sup>[1]</sup>, Csk<sup>[2]</sup>
```

Fluorescence stop-flow kinetics: release factor^[3], eIF5B^{[4],} EF-1B^[5], EF-Tu^[6]

FRET: AC^[7]

Specific Ligands:

Translation factor: eIF5B^[4], EF1B2^[5], EF-Tu^[6]

Release factor[^{3]}

Selected References:

[1]DnaB-DnaC protein complex with nucleotide cofactors. 1. Allosteric conformational transitions of the complex. *Biochemistry* **48**:6712.

[2] Shaffer *et al.* (2001) Nucleotide release and associated conformational changes regulate function in the COOH-terminal Src kinase, Csk. *Biochemistry* **40**:11149.

[3] Pisareva *et al.* (2006) Kinetic analysis of interaction of eukaryotic release factor 3 with guanine nucleotides. *J. Biol. Chem.* **281**:40224.

[4] Pisareva *et al.* (2007) Kinetic analysis of the interaction of guanine nucleotides with eukariotic translation initiation factor eIF5B.*Biochemistry* **46**:2622.

[5] Pittman *et al.* (2006) Mg2+ and key lysine modulate exchange activity of eukariotic translation elongation factoe 1B_. *J. Biol. Chem.***281**:19457.

[6] Wieden *et al.* (2002) Mechanism of elongation factor (EF)-Ts-catalyzed nucleotide exchange in EF-Tu. Contribution of contacts at the guanine base. *J. Biol. Chem.* **277**:6032.

[7] Goettle *et al.* (2007) Molecular analysis of the interaction of Bordetella pertussis adenylyl cyclase with fluorescent nucleotides. *Molecular Pharmacology* **72**:526.

Wehner et al. (2012) The guanine cap of human guanylate-binding protein 1 is responsible for dimerization and selfactivation of GTP hydrolysis. FEBS J. 279(2):203.

Oreb *et al.* (2011) Substrate binding disrupts dimerization and induces nucleotide exchange of the chloroplast GTPase Toc33.. *Biochem J.***436(2)**:313.

Koshiba *et al.* (2011) Structure-function analysis of the yeast mitochondrial Rho GTPase, Gem1p: implications for mitochondrial inheritance. *J. Biol. Chem.* **286(1)**:354.

Yang *et al.* (2010) A guaninine nucleotide exchange factor is a component of the meiotic spindle pole body in Schizosaccharomyces pombe. *Mol. Biol. Cell.* **21(7)**:1272.

Dahl *et al.* (2008) The Importance of P-loop and Domain Movements in EF-Tu for Guanine Nucleotide Exchange. *J. Biol. Chem.* **281(30)**:21139.

Ugolev *et al.* (2008) Dissociation of Rac1(GDP). RhoGDI Complexes by the Cooperative Action of Anionic Liposomes Containing Phosphatidylinositol 3,4,5-Trisphosphate, Rac Guanine Nucleotide Exchange Factor, and GTP. *J. Biol. Chem.* **283 (32)**:22257.

Ugolev *et al.* (2008) Liposomes Comprising Anionic but Not Neutral Phospholipids Cause Dissociation of Rac(1 or 2).RhoGDI Complexes and Support Amphiphile-independent NADPH Oxidase Activation by Such Complexes. *J. Biol. Chem.* **281 (28)**:19204.

Bos *et al.* (2007) Structure of the cyclic-AMP-responsive exchange factor Epac2 in its auto-inhibited state. *Nature* **439**:625.

Noel et al. (2007) Turning off the Ras switch with the flick of a finger. Nature Structural & Molecular Biology 4:677.

Pestova *et al.* (2006) Kinetic Analysis of Interaction of Eukaryotic Release Factor 3 with Guanine Nucleotides. *J. Biol. Chem.* **281** (52):40224.

Pick *et al.* (2006) Liposomes Comprising Anionic but not neutral Phospholipids cause dissociation of Rac(1 or 2)RhoGDI Complexes and support Amphiphile-independent NADPH oxidase Activation by Such Complexes . *J. Biol. Chem.* **281 (28)**:19204.

Wintermeyer *et al.* (2006) Role and timing of GTP binding and hydrolysis during EF-G-dependent RNA translocation on the ribosome. *PNAS* **103** (**37**):13670.

Thanbichler *et al.* (2000) Kinetics of the interaction of translation factor SelB from *Escherichia coli* with guanosine nucleotides and selenocysteine insertion sequence RNA. *J. Biol. Chem.* **275**:20458.

Ahmadian *et al.* (1997) Individual rate constants for the interaction of Ras proteins with GTPase-activating proteins determined by fluorescence spectroscopy. *Biochemistry* **36**:4535.

Remmers *et al.* (1996) Partial G protein activation by fluorescent guanine nucleotide analogs. Evidence for a triphosphatebound but inactive state. *J. Biol. Chem.* **271**:4791.

Rensland *et al.* (1995) Substrate and product structural requirements for binding of nucleotides to H-ras p21: the mechanism of discrimination between guanosine and adenosine nucleotides. *Biochemistry* **34**:593.

Richter *et al.* (1995) Interferon-induced MxA protein. GTP binding and GTP hydrolysis properties. *J. Biol. Chem.* **270**:13512.

Neal *et al.* (1990) Hydrolysis of GTP by p21NRAS, the NRAS protooncogene product, is accompanied by a conformational change in the wild-type protein: use of a single fluorescent probe at the catalytic site. *Proc. Natl. Acad. Sci. USA* **87**:3562.

Storage Conditions: store at -20 °C, protect from light

Precautions and Disclaimer:

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

©2016 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.