



QUALITY CONTROL RANGES

HUMAN FGF-21 ELISA CATALOG # EZHFGF21-19K

HUMAN FGF-21 STANDARD, LYOPHILIZED CATALOG# E8019-K

STANDARD	HUMAN FGF-21 VALUE
Lot Number. 8019-15K	2000 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6019-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	FGF21-115	66 – 137	pg/mL
QC 2	FGF21-215	511 – 1062	pg/mL