

1.09877.0010

Microbiology

OGYE Selective Supplement

Contents: 10 vials of lyophilisate for the preparation of 5 l
OGYE-Agar

Composition (per vial)

Oxytetracycline in a buffered base 0.05 g

Preparation of OGYE-Agar

1. Suspend 18.5 g OGYE-Agar (Base) in 500 ml of demin. water and autoclave (15 min. at 121 °C).
2. Let the medium cool to 45–50 °C.
3. Dissolve the supplement by adding 10 ml sterile distilled water and add the solution to the base medium. Mix gently and pour of plates.

Ordering information

1.05978.0500 OGYE-Agar (Base) 500 g

Version 2009-08-06

Merck KGaA, 64271 Darmstadt, Germany
Tel. +49 (0)61 51 72-24 40, www.merck-chemicals.com

