

53436 Perfringens S.F.P. Selective Supplement (S.F.P. Selective Supplement; Shahidi Ferguson Perfringens Selective Supplement)

An antimicrobial supplement recommended for the selective isolation of *Clostridium perfringens*.

Composition

(per vial, sufficient for 500 ml medium):

Kanamycin sulphate 6 mg

Polymyxin B 15000 IU

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 475 ml of sterile, molten Perfringens Agar Base (Prod No 39727), alongwith 25 ml of sterile Egg-Yolk Emulsion (Prod No 17148). Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

