



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP TGF β -1 Single Plex Magnetic Bead Panel Kit Catalog No. TGFBMAG-64K-01

Control Catalog No. LTGF-6064 Lot M664-113 and M664-213

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
TGF β 1	Control I	M664 - 113	161 - 335	pg/mL
	Control II	M664 - 213	556 - 1155	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

LTGF-6064: M664-113.213 (30 APRIL 2025) • SigmaAldrich.com/support • SigmaAldrich.com

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.