


NovAseptic[®] Valve Pneumatic Stainless Steel Actuator, Multiple Acting

Catalogue Number NA##/140-P0

Pneumatic Stainless Steel Actuator, Multiple Acting Configuration



Nominal Dimensions in mm (in.)

NA##	Size	A	B	C	D*	E	F	Compatible Valve Body	
								NA##	NU###
NA25	1"	133.0 (5.236)	73.0 (2.874)	84.0 (3.307)	TR75x2	M6	30.0 (1.181)	NA25	NU100
NA38	1½"	143.2 (5.638)	83.0 (3.268)	94.5 (3.720)	TR85x2	M6	33.0 (1.299)	NA38	NU150
NA51	2"	185.0 (7.283)	100.0 (3.937)	120.0 (4.724)	TR109x2	M6	33.0 (1.299)	NA51	NU200

* Note: Non-standard thread, Actuator specific.

Specifications

Net Weight (approximate)

Actuator	NA25	NA38	NA51
Weight kg (lb)	2.7 (5.95)	3.6 (7.94)	6.5 (14.3)

Material

	Cover Stainless Steel in Compliance with	
Material Code	EN 1.4305	ASTM® 303

Design Pressure

		NA25	NA38	NA51
Max	bar(g)	6.0	6.0	6.0
	(psi(g))	(87)	(87)	(87)
Min	bar(g)	-1.0	-1.0	-1.0
	(psi(g))	(-14)	(-14)	(-14)

Air Feed Pressure

		NA25	NA38	NA51
Max	bar(g)	7.0	7.0	7.0
	(psi(g))	(101)	(101)	(101)
Min	bar(g)	5.5	5.5	5.5
	(psi(g))	(80)	(80)	(80)
Pneumatic Fast Fitting mm (in.)		6 (0.24)	6 (0.24)	6 (0.24)

Specifications

Surface Roughness:	External surface Ra ≤ 1.6µm (63µin)
Environment Temperature:	0 to 70 °C (32 to 158 °F)
Note!	The applied valve body and diaphragm may have different design temperature and/or pressure limits. The weakest component in the assembled product determines the maximum design temperature and pressure limits.
Labeling:	Each Actuator is labeled for full LOT traceability according to Millipore QA routines.
Packaging:	The Actuator is packaged in a closed box.
Quality Control:	Millipore Quality Assurance System guarantees the control and traceability of the product.
Options:	For non-standard NovAseptic Actuator Options, please contact Millipore for further information.

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office or visit us at www.millipore.com.
In the U.S. and Canada, call toll-free 1-800-MILLIPORE (1-800-645-5476). In the U.S., Canada and Puerto Rico, fax orders to 1-800-MILLIFX (1-800-645-5439).

MILLIPORE

Millipore and NovAseptic are registered trademarks of Millipore Corporation or an affiliated company.

ASTM is a registered trademark of the American Society for Testing and Materials.

Specifications subject to change without notice.

Lit. No. SP1117EN00 Rev. 1, 09/2006. Replaces DS-061

© 2006 Millipore Corporation. All rights reserved.