


NovAseptic® Valve Tightening Tool for Locking Ring on NovAseptic Valve

Catalogue Number NA##/90

Tightening Tool for Locking Ring on NovAseptic Valve Configuration



Nominal Dimensions in mm (in.)

| NA##/90 | Compatible Valve Body NA | Compatible Valve Body NU |
|-----------|--------------------------|--------------------------|
| NA12/90 | NA12 | NU050 |
| NA18/90 | NA18 | NU075 |
| NA25/90 | NA25 | NU100 |
| NA38/90 | NA38 | NU150 |
| NA51/90 | NA51 | NU200 |
| NA76/90** | NA76 | NU300 |

* Socket for torquewrench compatible for NA12/90 through NA51/90, not included.
 ** The geometry of the NA76/90 has a different shape from the picture above, the tool is shaped as a standard hook wrench tool.

Specifications

Net Weight (approximate)

| Tightening Tool | NA12/90 | NA18/90 | NA25/90 | NA38/90 | NA51/90 | NA76/90 |
|-----------------|------------|------------|------------|------------|------------|-----------|
| Weight kg (lb) | 0.1 (0.22) | 0.1 (0.22) | 0.1 (0.22) | 0.1 (0.22) | 0.2 (0.44) | 0.8 (1.8) |

Material

| Tightening Tool | NA12/90 | NA18/90 | NA25/90 | NA38/90 | NA51/90 | NA76/90 |
|-----------------|-----------|---------|---------|---------|---------|---------|
| Material Code | Aluminium | | | | | Steel |


Surface Treatment: Anodized aluminium except on NA76/90.

Packaging: The Tightening Tool is packaged in a closed box.

Quality Control: Millipore Quality Assurance System guarantees the control at all stages of the manufacturing.

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office or visit us at www.millipore.com. In the U.S. and Canada, call toll-free 1-800-MILLIPORE (1-800-645-5476). In the U.S., Canada and Puerto Rico, fax orders to 1-800-MILLIFX (1-800-645-5439).



Socket for torquewrench installed on Tightening Tool for NovAseptic Valve, not included.



Millipore and NovAseptic are registered trademarks of Millipore Corporation or an affiliated company.

Specifications subject to change without notice.

Lit. No. SP1161EN00 Rev. -, 09/2006. Replaces DS-242

© 2006 Millipore Corporation. All rights reserved.