



ESGRO COMPLETE™ GELATIN SOLUTION

CATALOG NUMBER: SF008

QUANTITY: 500 mL

LOT NUMBER:

DESCRIPTION: ESGRO Complete Gelatin Solution is qualified for use with mouse embryonic stem cells cultured in serum-free conditions with ESGRO Complete PLUS Clonal Grade Medium (Cat. No. SF001-100P or SF001-500P). Product is ready to use for plating.

APPLICATIONS: Plate coating

PROTOCOL: **Gelatin coating protocol – Serum Free Media**

1. Warm 0.1% Gelatin Solution to room temperature before use.
2. Add enough Gelatin Solution to adequately cover plasticware surface.
3. Leave with lid on in laminar flow hood for 30 minutes.
4. Aspirate Gelatin Solution.
5. Add media and cell suspension to plasticware.

FORMAT: 0.1% gelatin solution sourced from porcine.

PRESENTATION: Liquid in sterile water.

STORAGE/HANDLING: Store at ambient temperature. Refer to lot expiration date on label.

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2010: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.