

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

## **ProductInformation**

## **BORIC ACID TABLETS**

Product No. B2645

## **Product Description**

Boric Acid is typically used as a preservative/stabilizer for urine specimens. Boric Acid in the tableted format also provides the advantage of a neat and convenient source of boric acid for general laboratory use.

Appearance White, round tablet, 5/8" in diameter Tablet Weight 1.0 gram (range 0.950 - 1.050 g)
Solubility\* Clear, colorless solution
Dissolution Time\* Not more than 30 minutes pH 5.0 (range 4.8 - 5.2)
Assay 100% (range 95-105%)

(NaOH Titration)

\*One tablet dispersed in 100 ml of deionized water.

Packaging 100 or 1000 tablets per bottle.

Custom packaging and bulk purchase information available

upon request.

Storage Store tablets at room temperature.

11/98