

shelf life information data sheet

Product Family: Polygard® CN

shelf life information

Shelf Life: 60 Months

Expiration Date: N/A

Storage Conditions: Ambient temperature

Applicable Catalog Number(s): See Attached list

quality attributes studied

Visual Examination

Water Flow Rate

Gravimetric Extractables

Not all products listed have undergone stability testing to determine the products' shelf lives. The design of the stability schedule is such that only select product(s) within a product family are tested. The design assumes that the products selected are representative of the family as they utilize the same materials of construction, production processes, and packaging components as the products within the family.

The product(s) above does not have a published expiration date or shelf life claim. The shelf life stated above is a guideline based upon available product stability data which is constantly updated as new data is made available. It is not a guarantee or a claim designating the time prior to which a lot of the product will remain within the shelf life specification if stored under the defined conditions, and after which it must not be put into use. It is not an expiration date or a use-by-date. Further, as this is not a product claim, updates to product shelf lives do not require customer notification.

For firms that must establish expiration dates or use-by-dates for this product(s) utilized in a cGMP manufacturing environment, it is recommended that the firms assign a default use-by-date based on their internal procedures/practices.

A handwritten signature in black ink, appearing to read "T. Brackett".

01-May-19

Thomas Brackett
Supplier Quality Engineering Manager
Third Party Finished Goods Quality

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.



CN0301E06	CN0601E06	CN1202E06	CN12M4T06	CN2501V06	CN3HM4E06
CN0301E06	CN0601E06	CN1202E06	CN12M4T06	CN2501V06	CN5001E06
CN0301T06	CN0601S06	CN1203S06	CN1H01E06	CN2502T06	CN5001T06
CN0301T06	CN0601T06	CN1203S06	CN1H01E06	CN2502T06	CN5001T06
CN0302E06	CN0602E06	CN1203V06	CN1H01V06	CN2571E06	CN5001V06
CN0302E06	CN0602E06	CN1203V06	CN1H01V06	CN2571S06	CN5001V06
CN0303E06	CN0603E06	CN1251S06	CN1H02E06	CN2571S06	CN5002E06
CN0303E06	CN0603E06	CN1251S06	CN1H02E06	CN2572E06	CN5002E06
CN0351E06	CN0603T06	CN1271E06	CN1H71E06	CN2572S06	CN5003E06
CN0351E06	CN0603T06	CN1271S06	CN1H71T06	CN2572S06	CN5003E06
CN0371E06	CN0671E06	CN1271S06	CN1H71V06	CN2573E06	CN5051E06
CN0371E06	CN0671E06	CN1271T06	CN1H71V06	CN2573S06	CN5051E06
CN0371S06	CN0671V06	CN1271T06	CN1H72E06	CN2573S06	CN5053E06
CN0371S06	CN0671V06	CN1271V06	CN1H72E06	CN25F1E06	CN5071E06
CN0372E06	CN0672E06	CN1271V06	CN1H72V06	CN25F2V06	CN5071S06
CN0372S06	CN0673E06	CN1272E06	CN1H72V06	CN25F2V06	CN5071T06
CN0372S06	CN0673E06	CN1272T06	CN1H73E06	CN25F3E06	CN5071T06
CN0373E06	CN0673S06	CN1272T06	CN1H73E06	CN25F3E06	CN5072E06
CN0373E06	CN0673S06	CN1272V06	CN1H73V06	CN25M4S06	CN5072E06
CN0373S06	CN0673V06	CN1272V06	CN1H73V06	CN25M4S06	CN5073E06
CN0373S06	CN0673V06	CN1273E06	CN1HF1E06	CN25M4V06	CN5073S06
CN03F1E06	CN06F1E06	CN1273S06	CN1HF1E06	CN3H01E06	CN5073V06
CN03F1E06	CN06F1E06	CN1273S06	CN1HF1S06	CN3H01E06	CN5073V06
CN03F2T06	CN06F3V06	CN1273V06	CN1HF1S06	CN3H02E06	CN5074T06
CN03F2T06	CN06F3V06	CN1273V06	CN1HF2E06	CN3H02E06	CN5074T06
CN03F2V06	CN06M4E06	CN1274V06	CN1HF2E06	CN3H71E06	CN5074V06
CN03F2V06	CN06M4E06	CN1274V06	CN1HM4E06	CN3H71E06	CN5074V06
CN03F3V06	CN06M4T06	CN12F1E06	CN2501E06	CN3H72E06	CN50F1E06
CN03F3V06	CN06M4T06	CN12F1S06	CN2501E06	CN3H73E06	CN50F1E06
CN03M0410	CN1201E06	CN12F1S06	CN2501S06	CN3H73E06	CN50F2E06
CN03M0410	CN1201E06	CN12F2E06	CN2501S06	CN3HF1V06	CN50F2E06
CN03M4E06	CN1201S06	CN12F2E06	CN2501T06	CN3HF1V06	CN50F3E06
CN03M4S06	CN1201T06	CN12M4E06	CN2501T06	CN3HM4E06	CN50M4E06
CN03M4S06	CN1201T06	CN12M4E06			