

eMcollect

Cat. No. A04551

Use

Sterile sampling device adapted to sample carcasses and work surfaces in food industries.

Used to check the environment during production when low probability or absence of disinfectant.

Pathogens detected: Listeria, Salmonella, other ones.

Features

Sponge

- Premium quality cellulose
- Size 8 x 4 cm and 1.2 to 1.6 cm thick
- Blue plastic made stick, with release clamps

Diluent

- Buffered peptone water
- Calibrated dose: 10 mL

Green pre-cut sampling bag

- Easy to open, 20 x 32.5 cm
- Large writing area for the sample information
- Adhesive closure label
- **Green resealable label** to close the bag after sampling, with marking "during production"

Product Data

Cat. No.	Custom code	Packaging	Cardboard box size
A04551	3921140000	100 units per box: 10 zip bags of 10 units	40 x 40 x 30 cm

Sterility and traceability:

Sterilization	Sterilization certificate	Shelf life	Traceability **
Sterilization by gamma rays Calibrated dose minimum 20 kiloGrays	in each box	12 months	On each product/box

** Printed information: Cat. No., designation, batch number, diluent code, production, and expiry date



Before



After





Cat. No. A04551

Sponge stick, BPW
in bag with closure label

Composition (ingredients in grams per liter of purified water):

- Casein peptone: 10.00
- Sodium chloride: 5.00
- Disodium phosphate, 12 H₂O: 9.00*
- Monopotassium phosphate: 1.50

* Corresponding to 3.56 g disodium phosphate anhydrous.



Green
resealable label
on the sampling
bag, for sampling
during production



Find the eMcollect brochure and further information on
our website **www.SigmaAldrich.com/HM** or scan this
code to get in touch with your local specialist

To place an order or receive technical assistance

Order/Customer Service: **SigmaAldrich.com/order**

Technical Service: **SigmaAldrich.com/techservice**

Safety-related Information: **SigmaAldrich.com/safetycenter**

SigmaAldrich.com

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
63843 | Rev. 09/2025



ISO 18593 standard
ISO 17604 standard