

Don't use printout for evaluation!!!




114411
MQuant®
Sulfate Test
Sulfat-Test
Test Sulfates
Test Sulfatos




Color card
Farbkarte
Carte colorimétrique
Tarjeta colorimétrica




1.



2.



3.



7.69144.1102-XXXXXX 06/2020

Merck KGaA, 64271 Darmstadt, Germany,
Tel. +49(0)6151 72-2440
www.sigmaldrich.com/mobile-laboratory

EMD Millipore Corporation, 400 Summit Drive
Burlington MA 01803, USA, Tel. +1-978-715-4321

Sigma-Aldrich Canada Co. or Millipore (Canada) Ltd.
2149 Winston Park, Dr. Oakville, Ontario, L6H 6J8
Phone: +1 800-565-1400

0

25

50

80

110

140

200

300

mg/l
 SO_4^{2-}

