



SIGMA QUALITY CONTROL TEST PROCEDURE

Product Information

Enzymatic Assay of ACYLASE I (EC 3.5.1.14) Sigma Prod. No. A-2156

PRINCIPLE:

N-Acetyl-L-Methionine + H₂O $\xrightarrow{\text{Acylase I}}$ L-Methionine + Acetic Acid

CONDITIONS: T = 37°C, pH = 8.0, A_{570nm}, Light path = 1 cm

METHOD: Colorimetric

REAGENTS:

- A. 100 mM N-Acetyl-L-Methionine Solution, pH 8.0 at 37°C (NAMET)
(Prepare 25 ml by dissolving 478 mg of N-Acetyl-L-Methionine, Sigma Prod. No. A-3258 in 10 ml of deionized water and 2 ml of 1 M NaOH. Adjust to pH 8.0 at 25°C with 10 M NaOH and then bring up to a volume of 25 ml with deionized water.)
- B. 200 mM Citrate Buffer, pH 5.0 at 37°C
(Prepare 200 ml in deionized water using Citric Acid, Free Acid, Anhydrous, Sigma Prod. No. C-0759. Adjust to pH 5.0 at 37°C with 1 M NaOH.)
- C. 1.6% (w/v) Stannous Chloride Solution (SnCl₂)
(Prepare 12 ml in Reagent B using Stannous Chloride, Anhydrous, Sigma Prod. No. S-2752.)
- D. Ethylene Glycol Monomethyl Ether
(Use Ethylene Glycol Monomethyl Ether, Sigma Prod. No. E-5378.)
- E. 2% (w/v) Ninhydrin Color Reagent (NCR)
(Prepare 50 ml by dissolving 1 g of Ninhydrin, Sigma Prod. No. N-4876, in 25 ml of Reagent D. Then add 25 ml of Reagent B. Store in an amber colored bottle.)
- F. 50% (v/v) 1-Propanol Solution
(Prepare 100 ml in deionized water using 1-Propanol, Sigma Stock No. 29,328-8.)

Enzymatic Assay of ACYLASE I
(EC 3.5.1.14)
Sigma Prod. No. A-2156

REAGENTS: (continued)

- G. 0.5 mM Cobalt Chloride Solution (CoCl_2)
(Prepare 20 ml in deionized water using Cobalt Chloride, Hexahydrate, Sigma Prod. No. C-2644.)
- H. 100 mM Barbitol Buffer, pH 8.0 at 37°C
(Prepare 25 ml in deionized water using Barbitol, Sodium Salt, Sigma Prod. No. B-0500. Adjust to pH 8.0 at 37°C with 1 M HCl.)
- I. 0.8 mM L-Methionine Standard Solution (Std Soln)
(Prepare 20 ml in deionized water using L-Methionine, Sigma Prod. No. M-9625.)
- J. Acylase I Enzyme Solution
(Immediately before use, prepare a solution containing 0.1 unit/ml of Acylase I in cold deionized water.)

PROCEDURE:

Step 1:

Pipette (in milliliters) the following reagents into a suitable container:

	<u>Test</u>
Reagent J (Enzyme Solution)	1.00
Reagent H (Barbitol Buffer)	2.00
Reagent G (CoCl_2)	1.00

Mix by swirling and equilibrate to 37°C. Then add:

Reagent A (NAMET)	1.00
-------------------	------

Immediately mix by swirling and incubate at 37°C for exactly 30 minutes. Remove 1 ml from the Test and place into a glass-stoppered test tube. Heat in a boiling water bath for 3 minutes.

Enzymatic Assay of ACYLASE I
(EC 3.5.1.14)
Sigma Prod. No. A-2156

PROCEDURE:

Step 2:

Pipette (in milliliters) the following reagents into a suitable container.

	<u>Blank</u>
Reagent J (Enzyme Solution)	1.00
Reagent H (Barbital Buffer)	2.00
Reagent G (CoCl ₂)	1.00

Mix by swirling and equilibrate to 37°C. Then add:

Reagent A (NAMET)	1.00
-------------------	------

Mix by swirling and immediately remove 1 ml from the Blank and place into a glass-stoppered test tube. Heat in a boiling water bath for 3 minutes.

Cool the Test and Blank solutions which have been boiled by using running water. Then add:

	<u>Test</u>	<u>Test Blank</u>	<u>Std 1</u>	<u>Std 2</u>	<u>Std 3</u>	<u>Std 4</u>	<u>Std 5</u>	<u>Std Blank</u>
Test Solution	1.00	----	----	----	----	----	----	----
Blank Solution	----	1.00	----	----	----	----	----	----
Reagent I (Std Soln)	----	----	0.10	0.20	0.40	0.60	1.00	----
Deionized Water	----	----	0.90	0.80	0.60	0.40	----	1.00
Reagent E (NCR)	2.00	2.00	2.00	2.00	2.00	2.00	2.00	2.00
Reagent C (SnCl ₂)	0.10	0.10	0.10	0.10	0.10	0.10	0.10	0.10

Mix by swirling and place the vials in a boiling water bath for 20 minutes. Remove the vials and allow to cool to room temperature. Add 10 ml of Reagent F (1-Propanol) to each vial. Mix well and transfer the vial contents to suitable cuvettes. Determine the absorbance at 570 nm for each of the vials using a suitable spectrophotometer.

Enzymatic Assay of ACYLASE I
(EC 3.5.1.14)
Sigma Prod. No. A-2156

CALCULATIONS:

Standard Curve:

$$\Delta A_{570\text{nm}} \text{ Standard} = A_{570\text{nm}} \text{ Standard} - A_{570\text{nm}} \text{ Standard Blank}$$

Prepare a standard curve by plotting the $\Delta A_{570\text{nm}}$ of the L-Methionine Standard Solution versus micromoles of L-Methionine.

Sample Determination:

$$\Delta A_{570\text{nm}} \text{ Sample} = A_{570\text{nm}} \text{ Test} - A_{570\text{nm}} \text{ Sample Blank}$$

Determine the μmoles of L-Methionine liberated using the Standard curve.

$$\text{Units/ml enzyme} = \frac{(\mu\text{moles of L-Methionine liberated})(5)(\text{df})}{(1)(1)}$$

5 = Total volume (in milliliters) of assay

df = Dilution factor

1 = Volume (in milliliter) of enzyme used

1 = Volume (in milliliter) of sample used in Colorimetric Determination

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

UNIT DEFINITION:

One unit will hydrolyze 1.0 μmole of N-acetyl-L-methionine per hour at pH 8.0 at 37°C.

FINAL ASSAY CONCENTRATION:

In a 5.00 ml reaction mix, the final concentrations are 40 mM barbital, 20 mM N-acetyl-L-methionine, 0.1 mM cobalt chloride, and 0.1 unit acylase I.

REFERENCE:

Moore, S. and Stein, W.H., (1948) *Journal of Biological Chemistry* **176**, 367-388

**Enzymatic Assay of ACYLASE I
(EC 3.5.1.14)
Sigma Prod. No. A-2156**

NOTES:

1. This assay is based on the cited reference.
2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.