

M-ENDO MEDIUM

Ampoules or vials of M-Endo Medium may vary in color from light pink to dark pink. Small specks of material may be found either adhering to the inside wall of the ampoule or vial or floating in the medium. The specks are from undissolved dye used in the raw medium formulation. The color variation and the dye specks do not affect either the growth recovery characteristics or the sheen production within the coliform colonies. The color of the non-coliform colonies is not affected.

NOTE: Add the medium directly to the absorbent pad beneath the membrane. Do not add the medium to the top of the membrane

CAUTION: This product is packaged in foil pouches to minimize exposure to light and oxygen, which can oxidize the medium and adversely affect its performance. For best results, store the product in sealed foil pouches at 4-8 C° until ready for use, and use up the contents of an opened pouch as quickly as possible.