

40952 Nickels and Leesment *ChromoSelect* Selective Supplement

An antibiotic supplement used for the selective isolation of citrate-fermenting lactic acid bacteria.

Composition (per vial, sufficient for 500 ml medium)
Vancomycin 100mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to sterile, molten, cooled (45-50°C) Nickels and Leesment *ChromoSelect* Medium (Cat. No. 68109). Mix well and pour into sterile Petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

