



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

LAMBDA PHAGE DNA, NON-METHYLATED

Product Number **D 3654**

Lot Number 071K9120

Storage Temperature: -20 °C

Product Summary

Source: Isolated from host strain *Escherichia coli* GM119. (rm⁻, dam⁻, dcm⁻)

Storage Buffer: 10 mM Tris-HCl,
pH 8.0, 1mM NaCl, 1 mM EDTA

Concentration: 400 µg/ml

A260/A280: 1.8

M.W.: 31.5 x 10⁶ Daltons

Base Pairs: 48,502

Unique Sites: Apa I, Nae I, Nar I, Nhe I, PaeR7 I,
SnaB I, Xba I, Xho I.

Precautions and Disclaimer

For laboratory use only. Not for drug, household or other uses.

Comments

Non-methylated DNA is completely digested by Bcl I, Cla I, Mbo I, Mbo II, Taq I or Xba I whereas methylated DNA (D 9768 or D 3779) is only partially cleaved by these enzymes.

Note: One A260 unit is equivalent to approximately 50 µg of DNA.

References

Sanger, F., et al., J. Mol. Biol. 162, 729 (1982).

9/01