



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone (800) 325-5832 (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## ProductInformation

### LAMBDA PHAGE DNA

Product No. **D3779**

Lot No. 69H9125

Storage: Below 0°C

### PRODUCT SUMMARY

**Source:** *Escherichia coli* W3110

**Buffer:** 10 mM Tris-HCl, pH 8.0, 1 mM EDTA

**Concentration:** 700 µg/ml or 14.0 A<sub>260</sub> units/ml

**Ratio:** A<sub>260</sub>/A<sub>280</sub>: 1.8

**M.W.:** 31.5 x 10<sup>6</sup> Daltons

**BASE PAIRS:** 48,502

**Unique Sites:** Apa I, Nae I, Nar I, Nhe I, PaeR7 I, SnaB I, Xba I, Xho I

### COMMENTS

Strain λc1857, Sam 7

Note: One A<sub>260</sub> unit is equivalent to approximately 50 µg of DNA.

**Reference:**

Sanger, F., et al., J. Mol. Biol., **162**, 729 (1982).

7/99

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.

Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.