

Product Information

Anti-Caveolin-1

Developed in Rabbit
Affinity Isolated Antibody

Product Number **C 4490**

Product Description

Anti-Caveolin-1 is developed in rabbit using as immunogen a synthetic peptide corresponding to a region at the N-terminus of human caveolin-1 (amino acids 2-20), conjugated to KLH. This sequence is highly conserved in many species, e.g., identical in rat, mouse, bovine, and dog caveolin-1. This sequence is not found in caveolin-2 or caveolin-3. Anti-Caveolin-1 is affinity-purified using the immunizing peptide immobilized on agarose.

Anti-Caveolin-1 recognizes caveolin-1. Applications include the detection and localization of caveolin-1 by immunoblotting (22 kDa) and immunofluorescence. Staining of caveolin-1 in immunoblotting is specifically inhibited with the caveolin-1 immunizing peptide (human, amino acids 2-20).

Caveolae are cholesterol/sphingolipid-rich, flask-shaped microdomains of the inner side of the plasma membrane with a diameter of 50 to 100 nm.^{1,2} Caveolae may play an important role in numerous essential cellular functions, including signaling, transport, lipid metabolism, cellular growth control and tumor suppression.³⁻⁵ Caveolae are present in many cell types, and are most abundant in endothelial cells, fibroblasts, smooth muscle cells and adipocytes.^{3,4} Caveolin, a 20-24 kDa integral transmembrane protein, has been identified as a principal component of caveolae membranes.

Caveolin was first identified as a major v-Src substrate in Rous sarcoma virus-transformed cells.⁶ Caveolin (also termed VIP21) exists in several isoforms termed caveolin-1, caveolin-2, and caveolin-3. Caveolin-1 (20-22 kDa) can exist as two isoforms, caveolin-1 α and -1 β due to alternative splicing of the mRNA. Caveolin-1 and -2 have similar tissue distribution. Caveolin-2 colocalizes with caveolin-1 and forms a hetero-

oligomeric complex with caveolin-1 *in vivo*.⁷ In contrast, caveolin-3 is a distinct isoform that is restricted to smooth, skeletal, and cardiac muscle.⁸ It has been proposed that the caveolin family members function as scaffolding proteins to organize and concentrate specific lipids such as cholesterol and glycosylphosphatidylinositol (GPI) and lipid modified signaling molecules within caveolae membranes.^{5,9,10} Caveolin can simultaneously recognize GPI-linked proteins and interact directly with a number of caveolae-associated downstream signaling molecules, such as H-Ras, hetero-trimeric G-proteins, annexin-II, EGF receptor, protein kinase C, src-family tyrosine kinases, and nitric oxide synthase (NOS) isoforms.^{5,11,12} Caveolae and caveolin-1 are both down-regulated in response to activated oncogenes such as H-Ras and v-Abl.¹³ In PC12 cells, caveolin-1 is upregulated during NGF-induced differentiation.¹⁴ Knockout mice lacking the caveolin-1 gene, also lack caveolae. The absence of this organelle impairs nitric oxide and calcium signaling in the cardiovascular system causing aberrations in endothelium-dependent relaxation, contractility and maintenance of myogenic tone.^{15,16} In addition, the lungs of knockout animals display thickening of alveolar septa due to uncontrolled endothelial cell proliferation and fibrosis, resulting in severe physical limitations in the caveolin-1 knockout mice.

Reagent

Anti-Caveolin-1 is supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide.

Antibody Concentration: 1.5-2.0 mg/ml

Precautions and Disclaimer

Due to the sodium azide content, a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

By immunoblotting, a working antibody concentration of 1-2 µg/ml is recommended using a whole cell extract of the human endothelial ECV304 cell line and the mouse fibroblast NIH3T3 cell line.

By indirect immunofluorescence, a working antibody concentration of 4-8 µg/ml is recommended using the HUVEC cell line.

Note: In order to obtain the best results using various techniques and preparations, we recommend determining the optimal working concentration by titration.

References

1. Lisanti, M.P., et al., *Mol. Membr. Biol.*, **12**, 121-124 (1995).
2. Okamoto, T., et al., *J. Biol. Chem.*, **273**, 5419-5422 (1998).
3. Smart, E.J., et al., *Mol. Cell. Biol.*, **19**, 7289-7304 (1999).
4. Li, S., et al., *J. Biol. Chem.*, **271**, 29182-29190 (1996).
5. Anderson, R.G., et al., *Ann. Rev. Biochem.*, **67**, 199-225 (1998).
6. Glenney, J.R., and Soppet, D., *Proc. Natl. Acad. Sci. USA*, **89**, 10517-10521 (1992).
7. Scherer, P.E., et al., *J. Biol. Chem.*, **272**, 29337-29346 (1997).
8. Song, K.S., et al., *J. Biol. Chem.*, **271**, 15160-15165 (1996).
9. Sargiacomo, M., et al., *Proc. Natl. Acad. Sci. USA*, **92**, 9407-9411 (1995).
10. Li, S., et al., *J. Biol. Chem.*, **271**, 568-573 (1996).
11. Couet, J., et al., *J. Biol. Chem.*, **272**, 30429-30438 (1997).
12. Oka, N., et al., *J. Biol. Chem.*, **272**, 33416-33421 (1997).
13. Engelman, J.A., et al., *J. Biol. Chem.*, **272**, 16374-16381 (1997).
14. Galbiati, F., et al., *Proc. Natl. Acad. Sci. USA*, **95**, 10257-10262 (1998).
15. Drab, M., et al., *Science*, **293**, 2449-2452 (2001).
16. Razani, B., et al., *J. Biol. Chem.*, **276**, 38121-38138 (2001).

KAA/ER 04/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.