

Product Name: SecureSeal™

Catalogue No.: T329-5

This document replaces any previous version

1. Product Description:

- **Aluminum Sealing Foil:** Disposable, non-sterile peeling foil for Microplates and PCR plates
- **Compatible with:**
 - T329-9 *Simport* Amplate™ Roller
 - T110 BioBlock™ Deep well plate series
 - T323 and T324 PCR plate series

2. Packaging:

- **Case:** 100 sheets per package

3. Product Specifications:

- **Material:**
 - Sealing Foil: Aluminum film with Acrylic adhesive
- Certified DNase, RNase, DNA Free and PR Inhibitor Free
- DMSO resistant
- Temperature range: -80°C to +120°

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **RoHS:** Plastics and colorants are in conformity with RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free

5. Quality Assurance:

- No presence of contamination in aluminum
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Check closure when using biohazard material and / or chemical reagents

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*