

**Product Name:** CLIKLOK™

**MICROCENTRIFUGE TUBE WITH FLAT CAP**

**Catalogue No.:** T330-6

**Material:**

Polypropylene Tube with attached cap

**Description:**

- 0.6 ml tube with attached cap
- Non sterile
- Conical bottom
- Graduated
- Etched surface on side of tube for sample identification
- Caps are pre-attached and not screwed on
- Caps have a loop and a unique ClkLok™ sealing mechanism
- Available color: Natural  
Blue, Green, Yellow on request only. Minimum quantity apply.

**Packaging:**

**Case:** 1 bag of 500 units

**Standards and Conformity:**

**Temperature range:** -80°C to +121°C

**Autoclavable:** +121°C for up to 30 minutes

**RNase, DNase, DNA and Pyrogen free:** certified

**ISO 2859-1:** Sampling and inspection

**Centrifugation:** up to 20 000g

**USP:** Plastic conforms to USP Class VI

**CONEG / RoHS:** Plastics comply

**REACH:** Plastic comply

**LATEX:** Plastic is Latex free

**BSE/TSE:** Plastic is BSE / TSE Free

**Recommended Use:** Designed for General laboratory use. Disposables for single use only. Verify proper cap closure when using biohazard material and / or chemical reagents. Follow chemical resistance chart recommendations.

**Storage Conditions:** Store at room temperature in normal warehouse conditions. Avoid temperature variations and humidity. Protect from any possible contamination and any damage to the packaging.

**Disclaimer:** *The information contained in this document shall not constitute a representation or warranty of the suitability for any specific purpose. To the best of our knowledge, the information contained is reliable and true based on our internal and external testing as well as statements supplied to us by third party suppliers.*