

Product Name: CLIKLOK™

MICROCENTRIFUGE TUBE WITH FLAT CAP

Catalogue No.: T330-6

Material:

Polypropylene Tube with attached cap

Description:

- 0.6 ml tube with attached cap
- Non sterile
- Conical bottom
- Graduated
- Etched surface on side of tube for sample identification
- Caps are pre-attached and not screwed on
- Caps have a loop and a unique ClikLok™ sealing mechanism
- Available color: Natural
Blue, Green, Yellow on request only. Minimum quantity apply.

Packaging:

Case: 1 bag of 500 units

Standards and Conformity:

Temperature range: -80°C to +121°C

Autoclavable: +121°C for up to 30 minutes

RNase, DNase, DNA and Pyrogen free: certified

ISO 2859-1: Sampling and inspection

Centrifugation: up to 20 000g

USP: Plastic conforms to USP Class VI

CONEG / RoHS: Plastics comply

REACH: Plastic comply

LATEX: Plastic is Latex free

BSE/TSE: Plastic is BSE / TSE Free

Recommended Use: Designed for General laboratory use. Disposables for single use only. Verify proper cap closure when using biohazard material and / or chemical reagents. Follow chemical resistance chart recommendations.

Storage Conditions: Store at room temperature in normal warehouse conditions. Avoid temperature variations and humidity. Protect from any possible contamination and any damage to the packaging.

Disclaimer: *The information contained in this document shall not constitute a representation or warranty of the suitability for any specific purpose. To the best of our knowledge, the information contained is reliable and true based on our internal and external testing as well as statements supplied to us by third party suppliers.*