

5-Carboxy-4'-aminomethylfluorescein

Catalog number: 207


Unit size: 10 mg

Product Details

Storage Conditions	Freeze (<15 °C), Minimize light exposure
Expiration Date	12 months upon receiving

Chemical Properties

Appearance	Orange solid
Molecular Weight	441.82
Soluble In	DMSO
Chemical Structure	



Spectral Properties

Excitation Wavelength	498 nm
Emission Wavelength	525 nm

Applications

5-Carboxy-4'-aminomethylfluorescein is a key building block for developing Big Dye DNA sequencing terminators. It can also be incorporated into a peptide sequence to give fluoresceinated peptides.