

6-TET, SE
[6-Carboxy-2',4,7',7-tetrachlorofluorescein,
succinimidyl ester]

Catalog number: 211
Unit size: 5 mg


Product Details

Storage Conditions	Freeze (<-15 °C), Minimize light exposure
Expiration Date	12 month upon receiving

Chemical Properties

Appearance	Red solid
Molecular Weight	611.17
Soluble In	DMSO

Chemical Structure



Spectral Properties

Excitation Wavelength	521 nm
Emission Wavelength	542 nm

Applications

6-TET, SE is a popular amino-reactive fluorescent probe that is widely used in nucleic acid sequencing and related research. It is also often used for labeling oligos.