

Product Name: Cryovial®

Catalogue No.: T301-5

Material:

Polypropylene Tube and Cap with Silicone seal

Description:

- 5 ml tube with cap
- Internal Threaded
- Round Bottom
- Printed with graduations and writing area
- Assembled with cap and red O-ring seal
- Compatible with color coding cap insert T312

Packaging:

Bags: 100 units

Case: 1000 units

Standards and Conformity:

Gamma Sterilized: SAL of 10^{-3}

Radiation specified dose: between 6.5 and 13.5 KGy

Temperature range: -196°C to +121°C

Autoclavable: +121°C for up to 30 minutes

RNase, DNase, DNA and Pyrogen free: certified

Leak-proof: tested in vacuum at 71.3cm Hg

Liquid Nitrogen resistant: Gas phase of LN²

ISO 2859-1: Sampling and inspection

FDA: Plastics conform to FDA 21 CFR 177.1520

USP: Plastic conforms to USP Class VI

CONEG / RoHS: Plastics comply

REACH: Plastic comply

LATEX: Plastic is Latex free

BSE/TSE: Plastic is BSE / TSE Free

MEA: Product has passed MEA testing

Recommended Use: Designed for storing biological material at temperatures as low as -196°C
Not recommended in the liquid phase of LN². Disposables for single use only.

Storage Conditions: Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging which could compromise the product sterility.

Disclaimer: *The information contained in this document shall not constitute a representation or warranty of the suitability for any specific purpose. To the best of our knowledge, the information contained is reliable and true based on our internal and external testing as well as statements supplied to us by third party suppliers.*