

Product Name: Cryovial® Workstation Rack

Catalogue No.: T315 Series

This document replaces any previous version

1. Product Description:

- **Cryovial® Workstation Rack:** Non-sterile, with locking base that can hold up to 50 Simport Cryovials® locked in place for one-hand operation
- **Available Colors:**
 - Blue, Red and Lilac

2. Packaging:

- **Case:** 4 units per case

3. Product Specifications:

- Material: Polypropylene
- Temperature range: -196°C to +121°C
- Autoclavable at 121°C, for up to 30 minutes
- Each well position identified by alphanumeric index

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **LATEX:** Material is Latex free

5. Quality Assurance:

- Clear, no presence of contamination in plastic;
- Visual attributes;
- Dimension measurements;
- Volume measurements.

6. Traceability:

- **Lot No. Composition:** 8 to 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label;
 2. On label inserted inside the master case;
 3. On the inner bag.

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions;
- Protect from any damage to the packaging.

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*